


NuLAB wet chemistry nutrient analyser

Verena Dauben¹, Alexander Davidov¹ & Vincent Kelly²

¹ MacArtney Germany GmbH, ² Green Eyes Environmental LLC, MD, USA

In 2018
MBT - Meerestechnisches Büro Turla GmbH
became
MacArtney Germany GmbH, Kiel

www.macartney.com



Key facts:

- located in Kiel
- part of the MacArtney Group – with over 400 employees and offices or partner on all continents

Co-operation with Green Eyes Environmental LLC:

- Since 2001 co-operation, service and support for former EnviroTech and later Green Eyes nutrient analysers
- Since 2016 sales representation of Green Eyes products in Europe
- 2017 – setting up of facilities at MacArtney Germany to test and service Green Eyes nutrient analysers and prepare reagents
- Nutrient analyser training, installations and workshops

Green Eyes Environmental LLC:

- Based in Easton, Maryland, United States
- Founded in 2006 by Vincent Kelly, Chemical Oceanographer
- Further development into state of the art analysers

Three NuLAB versions for various applications

Basic NuLAB

- smaller size (1-2 channels)




NuLAB Plus

- Included touch screen controller, relays for pump and water2web data posting
- Up to 3 channels



NuLAB Submersible (buoy or shallow water monitoring station)

- Up to 10 m water depth
- Up to 4 channels






- **What NuLAB does:**

- Application of established wet chemical methods (US EPA) to a field chemical analyser (for Nitrate, Phosphate, Ammonia and Silicate)
- Precise volumes of sample
- Data is calibrated via an On-Board-Standard (OBS)
- Reagents connected to a rotary valve and mixed by a syringe pump
- Analysed in high precision colorimeters.
- Operates with 8 pre-defined macros that determine how analyses are carried out

NuLAB Plus Set-up




Controller

Inlet with filter

**Reagent bags
(gas tight)**

Detector



8-port rotary valve

Syringe

Specification of NuLAB

Standard Ranges (detection limit to linear range, micro M)

mg/L:	N+N	Nitrite	Phosphate	Ammonium	Silicate
<u>High Sensitivity Detectors (2 mm)</u>	0.003 to 0.70	0.002 to 0.5	0.006 to 0.8	0.004 to 0.3	0.008 to 1.7
<u>Low Sensitivity Detectors (10 mm)</u>	0.01 to 2.8	0.008 to 2.1	0.025 to 2.0	0.02 to 1.0	0.04 to 2.8

micro mol/L:	N+N	Nitrite	Phosphate	Ammonium	Silicate
<u>High Sensitivity Detectors</u>	0.2 to 50	0.15 to 35	0.2 to 25	0.3 to 20	0.3 to 60
<u>Low Sensitivity Detectors</u>	0.8 - 200	0.6 - 150	1.0 - 75	1.5 to 75	1.5 to 100

- Precision (one SD @ midrange of scale): Nitrate 3%, Nitrite 2%, Phosphate 3%, Ammonium 3%, Silicate 3%
- Expanded Ranges: Up to 5 mg/l through auto-dilution
- Accuracy: Based on the accuracy of the preserved on-board standard and sample replicate precision
- Analyses: Typically 1000 per channel, per deployment. Controlled by reagent payload and chemistry
- Analysis Time: N+N 13 min, Nitrite 9 min, Phosphate 14 min, Ammonium 17 min, Silicate 16 min
- Consumption: Sample ~ 2 ml, each reagent ~ 0.1 ml and DIW


Analysis definition NuLAB

Aim: Turning mass equivalents into colours

- Macros can be customized to specific requirements and to third party analysis protocols i.e.
 - Change of the sample and reagent volumes
 - Change of the mixing times and volumes
 - Change of the flushing
 - Change of temperature and heating duration etc,

```
1 # ASL NITRATE STANDARD V1
2 # POWER ON
3 Y1
4 D2
5 # TURN ON HEATERS
6 H3
7 # ALIGN
8 U1
9 G1
10 # PORT H (AIR)
11 U1
12 P8
13 # RETRACT 1.0ML
14 U2
15 -15080
16 # PORT A (DETECTOR)
17 U1
18 P1
19 # INSERT 1.0ML
20 U2
21 +15080
22 # PORT F (STANDARD)
23 U1
24 P6
25 # RETRACT 0.25ML
26 U2
27 -3770
28 # PORT A (DETECTOR)
29 U1
30 P1
31 # INSERT 0.25ML
```



Manual Mode



NuLAB Software

- NuLAB is operated via a software running on the controller
- Software includes two different modes ("Manual" & "Logging")
- Analyses are executed by macros
- Macros can be selected via the software interface
- Individual channels can be selected
- Deployment mode is defined

Logging Mode



Development of the NuLAB for FerryBox:

- The NuLAB was originally designed for fixed station monitoring with sample intervals of two hours or longer.
- In contrast, FerryBox application often desire higher sample frequencies to resolve sharp spatial gradients

The following adjustments were introduced:

Elevated reaction temperatures:

- By raising reaction temperatures, the influence of ambient and sample water temperature on color development will be minimized.
- This will stabilize the on-board standard (OBS) results used for sample calibration and reduce the necessary OBS analysis frequency.
- This will reduce analysis time
- Elevating reaction temperature required simple hardware and software changes

Maintain detector heaters and LEDs for extended periods:

- To reduce detector warm up times, the heaters and LEDs can be left ON between analyses. This will also improve instrument precision.
- These changes are limited to the analytical macros and an additional “detector warm-up” command from the controller.

Addition of a solenoid switching valve:

- Currently all channels other than nitrate (when also measuring nitrite) have deionized water (DIW) connected to the eight-port rotary valve that is used for flushing and reagent blanks.
- Users are now able to add an optional solenoid valve to the nitrate channel so that all channels will be equipped with DIW.
- The valve will switch between imidazole buffer and DIW and be controllable via terminal commands, the NuLAB controller and NuLAB macros.
- This option requires an additional solenoid valve and relays

New controller:

- New microprocessor based controller allow one serial port of a computer or datalogger to run and collect data from up to four Nu-LAB channels simultaneously.
- This is better option than present Linux controller for FerryBox or other heavily integrated applications with master computers.

Conclusion: Strengths of NuLAB

- NuLAB determines a reagent blank before each sample
- An OBS can be measured before each sample and NuLAB uses the most recent OBS for concentration calculation
- NuLAB macros are customizable
- NuLAB is easy to operate and easy to integrate into other systems (e.g. Ferry Box)
- Data can be transmitted via internet to web gateways (like MetOcean Gateway)
- Support, development (improvement), compact, manageable but also affordable

Thank you for your attention!

Contacts:

Dr Verena Dauben – MBT GmbH


✉  ved@macartney.com

Dr Alexander Davidov – MBT GmbH

✉  ad@macartney.com



Chemical Methods used:

- **Ortho-phosphate P-PO₄ – Molybdenum reaction**
 - L. Drummond and W. Maher: Determination of phosphorus in aqueous solution via formation of the phosphoantimonymolybdenum blue complex - Reexamination of optimum conditions for the analysis of phosphate. *Analytica Chimica Acta* 302 (1995) pp. 69 – 74.
 - J. Murphy and J. P. Riley: A modified single solution method for the determination of phosphate in natural waters, *Analytical Chimica Acta*, 27 (1962) p. 31
- **Nitrate (N-NO₃) plus Nitrite (N-NO₂) Analysis**
 - J. D. H. Strickland and T. R. Parsons: *A Practical Handbook of Seawater Analysis*. Ottawa: Fisheries Research Board of Canada, Bulletin 167, 2nd Ed., 1972. 293 pp.
 - Grasshoff, K: *Methods of Seawater Analysis*, Verlag Chemie, Weinheim and New York, 1976, pp.149 – 156
- **Ammonium (NH₄) Analysis**
 - J. D. H. Strickland and T. R. Parsons: *A Practical Handbook of Seawater Analysis*. Ottawa: Fisheries Research Board of Canada, Bulletin 167, 2nd Ed., 1972. 293 pp.
 - L Solorzano: Determination of ammonia in natural waters by the phenolhypochlorite method. *Limnol. Oceanogr.* Vol.14(5). 1969. pp. 799-801.
- **Silicate (SO₄) Analysis**
 - Determination of Dissolved Silicate in Estuarine and Coastal Waters by Gas Segmented Continuous Flow Colorimetric Analysis
 - Jia-Zhong Zhang, Cooperative Institute for Marine and Atmospheric Studies, Rosenstiel School of Marine and Atmospheric Science.
 - Atlantic Oceanographic and Meteorological Laboratory, National Oceanic and Atmospheric Administration, University of Miami, Miami, FL 33149
 - George A. Berberian, National Oceanic and Atmospheric Administration, Atlantic Oceanographic and Meteorological Laboratory, Ocean Chemistry Division, Miami, FL 33149




NuLAB installation at Helgoland


- Installation of a 2-channel NuLAB system (Nitrate & Phosphate) next to Ferry Box based at "Helgoland Einlaufbauwerk"
- Continuous test measurements since mid-July 2017
- Hourly samples of NO_3+NO_2 and PO_4 from Ferry Box sampling water



NuLAB Nitrate+Nitrite (mg/L)
Tiffin, OH USA



NuLAB Reactive Phosphorus (mg-P/L)
Tiffin, OH USA



NuLAB Nitrate concentration Helgoland

